

Overview Pediatric Trauma Modules



With the add-on trauma modules, you can cover a wide range of trauma conditions to enhance training and simulations.

Pediatric Trauma Modules can be added to pediatric manikins to create childhood traumatic injuries, signs of abuse and scenarios to enhance training and simulations.

Module set includes:

Zygomatic lacerations	<ul style="list-style-type: none"> • Head skin with zygomatic lacerations, nasal bleeding • Head skin fits over manikin's head
Chest skin with 1st, 2nd and 3rd degree burns	<ul style="list-style-type: none"> • Burns illustrated with blistering in progressive sequence • Chest skin replaces manikin's chest skin and attaches at the shoulders • ECG monitoring / defibrillation not possible when module is in use
Road rash arm	<ul style="list-style-type: none"> • Left arm with road rash and fractured radius • Arm articulates to manikin at the shoulder
Arm with 1st, 2nd and 3rd degree burns	<ul style="list-style-type: none"> • Burns illustrated with blistering in progressive sequence • Left arm with plain forearm, articulates to manikin at the shoulder
Fractured radius / abnormal healing	<ul style="list-style-type: none"> • Lower right arm with abnormal healing of fractured radius • Right arm with plain upper arm, articulates to manikin at the shoulder
Dog bite hand	<ul style="list-style-type: none"> • Right hand with dog bite • Lower right arm, articulates to manikin at the elbow
Electrical burn entry	<ul style="list-style-type: none"> • Right hand with electrical wire burn entry • Lower right arm, articulates to manikin at

	the elbow
Tibia fracture	<ul style="list-style-type: none"> • Exposed fractured tibia protruding from lower leg • Lower right leg articulates to manikin at the knee
Closed femur fracture	<ul style="list-style-type: none"> • Left leg with closed femur fracture sleeve • Sleeve fits over manikin's left thigh
Cigarette burns	<ul style="list-style-type: none"> • Right leg with cigarette burns • Sleeve fits over manikin's right thigh
Crushed toes	<ul style="list-style-type: none"> • Left foot with crushed toes • Left lower leg articulates to manikin at the knee
Electrical burn exit	<ul style="list-style-type: none"> • Right foot with electrical wire burn exit wound • Right lower leg articulates to manikin at the knee