

AIRWAY MANAGEMENT TRAINER



Laerdal® Airway Management Trainer

Perfecting Procedures

Realistic practice is the key to developing proficiency in airway management skills. The Laerdal Airway Management Trainer realistically simulates an adult airway, and the lifelike upper torso and head can be used to demonstrate and practice intubation, ventilation, and suction techniques. In addition, it can be used to demonstrate upper airway bronchoscopy and bronchoscopy guided endotracheal intubation. Develop skills that are directly transferable to a clinical setting

Create an environment where learners can train with the equipment they use in real-life scenarios for suction and clearing, intubation, and bag-valve ventilation.

Bag-Valve Ventilation

- Provides visual inspection of lung expansion
- Provides auscultation of breath sounds
- Sellick Maneuver (Cricoid pressure)
- Laryngospasm

Intubation

- Tracheal (oral and nasal)
- Pharyngeal (oral and nasal)
- Esophageal
- Bronchial
- LMA insertion
- Endotracheal intubation
- Placement of supraglottic airways
- Correct tube placement
- Feedback for excessive laryngoscope pressure
- Bronchoscopic evaluation of tip position

Suction and Clearing

- Oral cavity
- Oro- or nasopharynx
- Oro- or nasotracheal, via endotracheal tube
- Gastric drainage